



January 13, 2009

RIT113 Required Configuration for Version 5.2:

1. Pentium-based or Centrino Duo computer (desktop or portable)
2. Adaptec SCSI Interface (for use with film digitizers)
(Adaptec 2930U for PCI, Adaptec 1460D PCMCIA, or Adaptec 1480B PCMCIA)
Or a USB2 interface **if using a Vidar Dosimetry Advantage Pro film scanner with USB2 capability**: Adaptec USB2Connect 4000 (full size PCI for desktop), or Adaptec USB2Connect (AUA-1420A) PCMCIA card for laptops.
3. Minimum RAM: 1 Gigabyte
4. Minimum hard disk space for installation: 1.1 Gigabyte.
5. CD-ROM Reader
6. A video graphics card and monitor capable of at least 1280x1024 resolution and capable of displaying at least 16 bit colors.
7. Operating System:
Windows 2000 Professional with Service Pack 4
or Windows XP Professional with Service Pack 2
8. Printer: Color or Black and White
(most facilities prefer color since many of the RIT113 outputs are in color)
9. Film Scanner:
 - Vidar VXR16 Dosimetry Pro with Firmware Version 13.0
 - Vidar Dosimetry Pro Advantage (Vidar part number 16321-004 or -005) with Firmware Version 46.3 and pass code suffix "BVVVNN"
 - Vidar Dosimetry Pro Advantage (Vidar part number 16321-007) (White LED) with Firmware Version 49.6, pass code suffix "BVVVNN" , serial number 341000 or greater.
 - Vidar Dosimetry Pro Advantage (Vidar part number 17500-001) (RED LED) with Firmware Version 49.7, pass code suffix "BVVVNN" , serial number 341000 or greater.
10. Floppy disk drive or thumb USB drive – only used if a license transfer is necessary.
11. Adobe® Reader – for viewing User's manual