



RIT®

TRACK & TREND WITH RIT® SOFTWARE

RIT TREND™ STATISTICAL DATABASE ANALYSIS

RITtrend™

The comprehensive statistical database solution that tracks and trends data for all of your department's measurements.

RIT TREND™ FEATURES

- Easily import your analysis results into the database with a single click.
- Set baseline measurements for statistical process control.
- Conveniently delete any data points that are not needed.
- Track historical data over time to detect any anomalies in your measurements.
- Integration with the RIT Tolerance Manager: Set fully customizable parameters to track and trend any pass/fail data in your set tolerances.
- Enable alerts via email, allowing you to immediately view significant changes in data or measurements that may exceed their set tolerances.



HOW TO USE RIT TREND™



Your data is automatically compared to your specifications with alerts for any discrepancies.

- ① Run the test analysis routine. (Using RIT or other software)
- ② With one click, add the analysis record to the RITtrend™ database or use the Multisource Data Manager to input data generated from non-RIT software.
- ③ If desired, create unlimited custom filtering specifications.
- ④ Use query builder selections to extract your data using customizable queries to sort and report multiple criteria, such as:
 - Analysis routine
 - Date
 - Equipment
 - User
 - Any predetermined custom fields, set by the user
- ⑤ Output Format Choices: View a visual representation of your data in a statistical process control chart, tabular format, or in customizable RITG142 and RITG135 reports.

REQUEST A PERSONAL DEMONSTRATION:

www2.radimage.com/demo

RITtrend™ is a trademark of Radiological Imaging Technology, Inc.